VascTM BAND hemostat

Instructions for Use

English/Instructions for Use	3
Français/Mode d'emploi	7

Vasc Band Hemostat

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.


CAUTION

The Vasc Band procedure should be performed by physicians or physician-directed allied health care professionals with adequate training in the use of the device.

DEVICE DESCRIPTION

The Vasc Band applies mechanical force at a percutaneous access site to achieve hemostasis. Each device consists of the following components:

- Compression band, consisting of:
 1. Adjustable retention strap
 2. Inflatable compression balloon
 3. Inflation valve
 4. Inflation tubing
- Inflation syringe (22ml)

The Vasc Band has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS – Vasc Band Hemostat

Model Number	3521	3524	3527	3529	3537
Band Size	21cm	24cm	27cm	29cm	37cm

INDICATIONS

The Vasc Band hemostat is a compression device to assist hemostasis of arterial, venous and hemodialysis percutaneous access sites.

CONTRAINDICATIONS

- Patients hypersensitive to the materials of compression device.
- Patients with infection or other serious skin diseases at the site of puncture.
- For radial artery catheterization only: patients with an abnormal Allen test or radial pulse, or insufficient dual arterial supply.
- Not intended for femoral artery compression.

WARNINGS

- Patients should not be left unattended while the Vasc Band is in use.
- Ensure correct alignment of the Vasc Band prior to use.
- Do not inject air into any port other than the air injection port of this device.
- Do not leave Vasc Band on for inappropriately long periods of time as tissue damage may occur.
- Arterial pulse distal to the compression site should be monitored to ensure the artery is not completely occluded as arterial damage and/or thrombosis could occur.

PRECAUTIONS

- Ensure correct placement of Vasc Band.
- The product is intended for single use only.
- Do not reuse or re-sterilize.
- Do not use if package is opened or damaged.
- Over inflation of balloon >18ml could cause balloon damage, compromising the performance of the Vasc Band.

ADVERSE EVENTS

Potential adverse events include allergic reaction, vessel occlusion, local venous thrombosis, hematoma formation, recurrent bleeding, peripheral nerve damage, and complex regional pain syndrome.

CLINICAL PROCEDURE

Carefully inspect the Vasc Band packaging and components for damage prior to use.

It is important to properly size the Vasc Band to the patient's extremity. There are five sizes available to accommodate different patients.

Application of the Vasc Band:

1. Using sterile technique, open the pouch and transfer the Vasc Band and syringe into the sterile field.
2. At the end of the catheterization procedure, withdraw the introducer sheath 2-3cm from the puncture site.

3. Align the center of the Vasc Band balloon 2-5mm proximal to the puncture site and wrap the Vasc Band straps around the patient's extremity. Secure the device to the patient's extremity using the hook and loop fastener. For optimal fit, secure the band around the extremity, allowing no room for slack, but do not overtighten the band.
4. Draw 18ml of air into the Vasc Band inflation syringe and connect the syringe to the inflation valve on the Vasc Band.
5. Slowly inject air into the Vasc Band balloon while simultaneously removing the sheath from the vessel. Once the sheath is completely removed, continue to inject air into the balloon until hemostasis is achieved. Note: The nominal air injection volume for the Vasc Band is 15ml and the maximum air injection volume is 18ml.
6. Slowly withdraw air from the Vasc Band balloon until there is oozing from the puncture site. Once oozing is observed, re-inject 2ml of air into the Vasc Band until complete hemostasis is achieved.
7. Record the patient's vital signs and assess perfusion by means of color, sensation, and temperature q5 minutes x 3. After the first 15 minutes, record vital signs and assess perfusion q15 minutes x 3, then q30 minutes until initiating device removal.
8. The Vasc Band should be left on the patient's extremity for two hours following catheterization unless directed otherwise by physician's orders or patient conditions.

Device Removal:

1. At the recommended device removal time, withdraw 2ml of air from the Vasc Band and observe the puncture site for bleeding.
2. If bleeding is present, inject enough air to restore hemostasis. Wait 30 minutes and then repeat step 1 of the removal protocol.
3. If no bleeding is present, continue to remove 2ml of air q15 minutes x 3 or until pressure is fully released.
4. When pressure is released and hemostasis confirmed, carefully remove the Vasc Band from the puncture site, being careful not to disrupt the clot. Apply dressing of choice.

The Vasc Band should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Note: If needed, a sizing aid is available from Vascular Solutions. To use the sizing aid:

1. With the colored side of the aid visible, place the blue START end of the aid at the puncture site.
2. Wrap the sizing aid snugly around the extremity, bringing the sizing line over the top of the measurement section.

3. The sizing line will indicate the optimal size Vasc Band to be used on the patient. If the sizing line is between two sizes, or if the indicated size is not available, use the larger size Vasc Band as indicated per the sizing aid.
4. Discard after use.

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the Vasc Band hemostat is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special or consequential damages arising from the use of the Vasc Band hemostat device. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.
See www.vasc.com/products/patents

Vasc Band is a trademark of Vascular Solutions, Inc.
See the International Symbols Glossary on page 11.

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Vasc Band est une marque de commerce de Vascular Solutions, Inc.